VULPES GOODS
BABYCARE

MATERNITY BAND

**MANUAL • ANLEITUNG
HANDLEIDING • MANUEL
MANUAL • MANUALE • MANUAL**

WWW.VULPESGOODS.COM

ENGLISH	4
The Device	5
User Guide	6
Characteristics and Application	10
DEUTSCH	11
Das Gerät	12
Benutzerhandbuch	13
Eigenschaften und Anwendung	17
NEDERLANDS	18
Het apparaat	19
Gebruiksaanwijzing	20
Kenmerken en Toepassingen	24
FRANCAIS	25
Le dispositif	26
Manuel de l'utilisateur	27
Caractéristiques et applications	31

ESPAÑOL

32

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ENGLISH

THE DEVICE

The maternity support belt has two sides: side A and side B. Side B is the inside and contains the printed text. This side should always be worn against the back.

Side A is the outside and is worn on the outside of the body. Both sides contain hook-and-loop fastening areas. The hook-and-loop parts on side B must be attached to the fastening areas on side A. When the inside (B) is put on first and the outside (A) is fastened over it, both sides form a firm and comfortable maternity support belt.

USER GUIDE

Step 1

Place the maternity belt low around the lower back and lower abdomen. Wrap the belt around the body and close it with the hook-and-loop fastener.

Step 2

Place the belly band as low as possible under the belly, without it rubbing against the legs while walking or sitting.

Step 3

Pull the belly band forward and attach it to the waist band. Fasten it securely with the hook-and-loop closure.

Step 4

The maternity support belt is now correctly attached and ready for use.

Usage instructions

- It is recommended to wear the belt over clothing. If you are comfortable with synthetic fabrics on the skin, you may also wear the belt under clothing. This is less visible and easier to match.
- The belt is best used while standing or walking. If you expect to sit for a long period, adjust the belt to the most comfortable position.
- When using the toilet, you may choose to detach the support band or leave it in place. If you choose not to detach it, lift the belly band above the belly and place it back after use.
- The fabric is soft and comfortable. Make sure the hooks of the fastening system do not attach directly to the fabric to prevent damage.
- If this happens accidentally, slight pilling or fuzzing is normal and will not affect the product's performance.

Washing instructions

- Fasten all hook-and-loop strips before washing to prevent them from attaching to other clothing and causing damage.
- Hand-wash the product in warm water with a mild detergent.
- Let it air dry. Do not iron or tumble-dry.
- If machine washing is necessary, use the laundry bag included with the product. Roll up the belt, make sure all hook-and-loop fasteners are closed and use a gentle washing program.
- Do not use harsh detergents or bleach.

CHARACTERISTICS AND APPLICATION

- The maternity belt is intended solely for support during pregnancy. Do not use this product for medical treatment or as a substitute for professional advice.
- Do not use the belt while sleeping, as this can create unnecessary pressure on the belly or lower back.
- Do not tighten the belt too much. Excessive pressure can cause discomfort or restrict circulation.
- Do not use the belt during intense sports activities or heavy physical effort.
- Do not wear the belt if you experience skin irritation, numbness or pain. Stop use immediately if discomfort occurs.
- The belt must not be used as a lifting aid, support strap or weight-bearing device.
- Avoid use in very warm conditions; excessive heat may cause discomfort.
- Do not allow small children to play with the belt due to choking or entanglement risks.
- The fastening system is strong; be careful that clothing does not get caught.
- Check regularly whether all parts are still securely attached and not worn out.

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